R	1 × 60 mL - Phosphate buffer pH 7 Phenol / Amino-4-antipyrine Cholesterol oxidase / Peroxydase
C	1 × 2 mL - Control solution

### Sample preparation

The concentration of Cholesterol in the sample used in the assay procedure has to be between 0.05 to 1 g/l.

### Precision

Under the described conditions, measurement accuracy is 5% on a control solution

### Assay procedure

Wavelength: 500nm / Optical path: 1cm / Temperature: 20-25 / 37°C  
Measurement: against water or air.

	Blank	Sample
Reagent	2,0 ml	2,0 ml
Mix and read DO	DO1	DO1
Water	100 µl	0
Sample	0	100 µl
Mix and read DO after 20 min	DO2	DO2

### Calculation

Determine the values for blank and samples:

$$\Delta DO = [ DO 2 - DO 1 ]_{\text{sample}} - [ DO 2 - DO 1 ]_{\text{blank}}$$

Cholesterol concentration is given by:

$$C = \frac{V \times MW}{\varepsilon \times l \times v \times 1000} \times \Delta DO \quad (\text{g/L})$$

In the assay procedure:

$$C = 1,1598 \times \Delta DO \quad [\text{g/l of Cholesterol in the sample}]$$

Dilution factor of the sample has to be considered in the calculation.

### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

### Assay control

Cholesterol reagents must be validated with the use of the include control.

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