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Sample preparation:

The concentration of citrate in the sample used in the assay procedure has to be between 0.04 to 0.4 g/l.

Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°
Measurement: against water or air.

	Blank	Sample
R 1	1 ml	1 ml
R 2	0.2 ml	0.2 ml
R 3	0.02 ml	0.02 ml
Water	1 ml	0.9 ml
Sample	0	0.1ml
Mix and read DO	DO 1	DO 1
R 4	0,02 ml	0,02 ml
Mix and read DO after 10 mn	DO 2	DO 2

Calculation

Determine the values for blank and samples:

$$\Delta DO = [DO 1 - DO 2]_{\text{échantillon}} - [DO 1 - DO 2]_{\text{blanc}}$$

Citrate concentration is given by:

$$C = \frac{V \times MW}{\epsilon \times l \times v \times 1000} \times \Delta DO \quad (g/l)$$

In the assay procedure:

$$C = 0,68302 \times \Delta DO \quad [g/l \text{ of citrate in the sample}]$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

The dissolved R4 is stable about 5 days, if stored at 2-8 °C.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.
Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 7.8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

Assay control

Citrate reagents must be validated with the use of the standard included in the kit. The standard is ready-to-use.

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