
# Ethanal kit

## Kit Ethanal

Cat. No. 038

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v038-21031

### Sample preparation:

The concentration of Ethanal in the sample used in the assay procedure has to be between 0.010 to 0.200 g/l.

Use clear, colorless and practically neutral liquid samples directly, or after dilution. Filter or centrifuge turbid solutions.

Ethanal is extremely volatile (boiling point = 21°C) :

- It is necessary to close tightly all containers with ethanol
- It is necessary to pipette solutions under the surface of the ethanal containing solutions

### Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

### Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°

Measurement: against water or air.

	Blank	Sample
R 1	1,0 ml	1,0 ml
R 2	0,2 ml	0,2 ml
Water	1,0 ml	0,9 ml
Sample	0	0,1 ml
Mix and read the DO	DO 1	DO 1
R 3	0,02 ml	0,02 ml
Mix and read DO after 10 min	DO 2	DO 2

### Calculation

Determine the values for samples:

$$\Delta\text{DO} = [\text{DO 2} - \text{DO 1}]_{\text{sample}} - [\text{DO 2} - \text{DO 1}]_{\text{blank}}$$

Ethanal concentration is given by:

$$C = \frac{V \times MW}{\varepsilon \times l \times v \times 1000} \times \Delta\text{DO Ethanal} \quad (\text{g/L})$$

In the assay procedure:

$$C = 0,15722 \times \Delta\text{DO} \quad [\text{g/l of Ethanal in the sample}]$$

Dilution factor of the sample has to be considered in the calculation.

### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8°C, contamination is avoided.

### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

### Assay control

Ethanal reagents must be validated with the use of the standard included in the kit.

Prepare one control at 0,200 g/L of ethanal.

Weight 29,0 mg of ethanal, and dissolve it in 100 ml water.

Solution is stable one day. Tightly close the flask.

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