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Preparation of the standard solution:

It is necessary to prepare a standard solution to dilute before use

Use one of the following products and to weigh in a precise way in a graduated flask of 100 ml the quantity indicated

Sodium sulphite, anhydrous, AR grade (Na ₂ SO ₃ , M=126.04 g/mol, approx 50.8% SO ₂)	60 mg
Sodium bisulphite, AR grade (Na ₂ S ₂ O ₅ , M=190.10 g/mol, approx 67.4% SO ₂)	45 mg
Potassium bisulphite, AR grade (K ₂ S ₂ O ₅ , M=222.33 g/mol, approx 57.6% SO ₂)	50 mg

That corresponds to 300 mg/L of total SO₂.

Dilute this solution to the 1/6: 1.5 ml + 7.5 ml of demineralised water. The total SO₂ concentration is then 50 mg/L.

Preparation of the sample:

The total SO₂ concentration in the sample used for the test has to be between 5 and 300 mg/L.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 405nm / Optical path: 1cm / Temperature: 20-37°
Measurement: against water or air.

	Blank	Standard	Sample
R 1	0,125 ml	0,125 ml	0,125 ml
Water	1,250 ml	1,235 ml	1,235 ml
Standard	0	0,015 ml	0
Sample	0	0	0,015 ml
Mix and read the DO	DO 1	DO 1	DO 1
R 2	0,125 ml	0,125 ml	0,125 ml
Mix and read the DO after 20 mn	DO 2	DO 2	DO 2

Calculation

Determine the following values for each cuvette
ΔDO = DO 2 – DO 1

Determine the absorbance differences then

For standard : $\Delta A_{\text{standard}} = \Delta DO_{\text{standard}} - \Delta DO_{\text{blank}}$
For sample : $\Delta A_{\text{sample}} = \Delta DO_{\text{sample}} - \Delta DO_{\text{blank}}$

Determine the total SO₂ concentration in the samples as described below

$$C(\text{mg/L of total SO}_2) = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{s tan dard}}} \times C_{\text{s tan dard}}$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry.
Storage conditions: 18 – 25 °C, darkness.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.
Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

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