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By mixing sample and chromogen in a pH 4.8 buffered solution, the reaction produces a measurable blue color within few minutes.

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The R3 reagent must be dissolved in 20 ml of R1 reagent. This solution is stable 3 week with 4°C.

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**Preparation of the standard solution:**

It is necessary to prepare a standard solution, included in the kit

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**Preparation of the sample:**

The iron concentration in the sample used for the test has to be between 0.1 and 15 mg/L.

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**Precision:**

Under the described conditions, measurement accuracy is 5% on a control solution

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**Assay procedure:**

Wavelength: 590nm / Optical path: 1cm / Temperature: 20-37°  
Measurement: against water or air.

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**Calculation:**

Determine the following values for each cuvette

$$\Delta DO = DO 2 - DO 1$$

Then, determine the absorbance differences

For standard :  $\Delta A_{\text{standard}} = \Delta DO_{\text{standard}} - \Delta DO_{\text{blank}}$

For sample :  $\Delta A_{\text{sample}} = \Delta DO_{\text{sample}} - \Delta DO_{\text{blank}}$

Finally, determine the Iron concentration in the samples as described below

$$C_{\text{iron}} (\text{mg} / \text{L}) = \left[ \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} \right] \times F$$

Where F = dilution factor of the sample (if diluted)

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**Storage instructions and reagent stability:**

The reagents are stable up to the indicated date of expiry.  
Storage conditions: 18 – 25 °C, darkness.

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**Warnings and precautions:**

Do not swallow. Avoid contact with the skin and mucous membranes.  
Take necessary precautions for the use of laboratory reagents.

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