R1	1 x 30 mL - Buffer pH 7,5 / NADP / ATP
R2	1 x 0,6 mL - Hexokinase / Glucose-6-phosphate dehydrogenase
RH	1 x 0,6 mL - Trehalase
CGlc	1 x 2 mL - Glucose Control solution Solution de contrôle de Glucose
CTre	1 x 1 g - Trehalose Control

Sample preparation:

The concentration of Trehalose in the sample used in the assay procedure has to be between 0.1 to 1 g/l.
If necessary, extract trehalose from sample with hot water (80°C)
(Ex.: Weigh 1.0 g of sample into a beaker, and add 40 ml of hot distilled water (~80°C). Mix 15 min. Transfer the solution into a 50 ml volumetric flask and adjust the volume)
Filter if necessary.

Removal of reducing sugars:

If sample contain glucose, it is necessary to remove it.
- Take 0.2 ml of the solution to be analysed
- Add 0.2 ml of a solution of sodium Borohydrate at 10 mg/ml in 50 mM sodium hydroxide. Mix vigorously and store at 40°C for 30 min
- Add 0.5 ml of a solution of acetic acid at 200 mM. Mix vigorously with a vortex mixer. *One effervescence should be observed.*

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°C
Measurement: against water or air.

	Blank	Sample
R 1	1,0 ml	1,0 ml
Water	2,0 ml	1,9 ml
Sample	0	0,1 ml
R 2	0,02 ml	0,02 ml
Mix and read the DO after 10 min	DO1	DO 1
R H	0,02 ml	0,02 ml
Mix and read DO after 20 min	DO 2	DO 2

Calculation

Determine the values for blank and samples:
 $\Delta DO = [DO 2 - DO 1]_{\text{sample}} - [DO 2 - DO 1]_{\text{blank}}$

Trehalose concentration is given by:

$$C = \frac{V \times MW}{\varepsilon \times l \times v \times 1000 \times 2} \times \Delta T_{\text{tréhalose}} \text{ (g/L)}$$

In the assay procedure:

$$C = 0,8283 \times \Delta DO \text{ [g/l of Trehalose in the sample]}$$

Dilution factor of the sample has to be considered in the calculation.
If reducing sugars had been removal, the result must be multiplied by the factor 0.9/0.2, by 4.5

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.
Take necessary precautions for the use of laboratory reagents.

Assay control

Trehalose reagents must be validated with the use of standard included in the kit. The Glucose standard is ready-to-use.
For Trehalose, weigh 500 mg by liter : it corresponds to a solution at 447 mg/L of Trehalose.

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