Enzymatic UV 340nm test

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Sample preparation:

The concentration of acetic acid in the sample used in the assay procedure has to be between 0.03 to 0.15 g/l.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°
Measurement: against water or air.

	Blank	Sample
R 1	1,0 ml	1,0 ml
R 2	0,2 ml	0,2 ml
Water	2,0 ml	1,9 ml
Sample	0	0,1 ml
Mix and read DO	DO1	DO 1
R 3	0,01 ml	0,01 ml
Mix and read DO after 3 min Start reaction with:	DO 2	DO 2
R 4	0,02 ml	0,02 ml
Mix and read DO after 10-15 min	DO 3	DO 3

Calculation

Determine the values for blank and samples:

$$\Delta A1 = [DO 2 - DO 1] \text{ et } \Delta A2 = [DO 3 - DO 1]$$

Then, use the following formula :

$$\Delta \text{Acetate} = \left[\Delta A2_{\text{sample}} - \frac{\Delta A1_{\text{sample}}^2}{\Delta A2_{\text{sample}}} \right] - \left[\Delta A2_{\text{blank}} - \frac{\Delta A1_{\text{blank}}^2}{\Delta A2_{\text{blank}}} \right]$$

Acetate concentration is given by:

$$C = \frac{V \times MW}{\varepsilon \times l \times v \times 1000} \times \Delta \text{acetate} \text{ (g/L)}$$

In the assay procedure:

$$C = 0,3079 \times \Delta \text{acetate} \text{ [g/l of acetic acid in the sample]}$$

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Once reconstituted, reagent R2 is stable 4 week

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.
Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

Assay control

Acetic acid reagents must be validated with the use of the standard included in the kit. The standard is ready-to-use.

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