R1	1 × 30 mL – Buffer pH10 / Glutamate 440mg
R2	1 × 6 mL – NAD 210mg
R3	1 × 0.6 mL – GPT 1100U
R4	1 × 0.6 mL – D-LDH 3800U
R5	1 × 0.6 mL – L-LDH 3800U
C	1 × 2 mL - Control solution Solution de contrôle

Sample preparation:

The concentration of D and L-Lactic acids in the sample used in the assay procedure has to be between 0.03 to 0.2 g/l.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°
 Measurement: against water or air.

	Blank	Sample
R1	1 mL	1 mL
R2	0.2 mL	0.2 mL
R3	0.02 mL	0.02 mL
Water	1 mL	0.9 mL
Sample	0	0.1 mL
Mix and read DO	DO1	DO1
R4	0.02 mL	0.02 mL
Mix and read DO after 45mn	DO2	DO2
R5	0.02 mL	0.02 mL
Mix and read DO after 45mn	DO3	DO3

Calculation

Determine the values for blank and samples:

For D-lactic acid:

$$\Delta DO_{D-lac} = [DO 2 - DO 1]_{\text{sample}} - [DO 2 - DO 1]_{\text{blank}}$$

For L-lactic acid:

$$\Delta DO_{L-lac} = [DO 3 - DO 2]_{\text{sample}} - [DO 3 - DO 2]_{\text{blank}}$$

Lactate concentration is given by:

$$C = \frac{V \times MW}{\epsilon \times l \times v \times 1000} \times \Delta DO \quad (\text{g/L})$$

In the assay procedure:

$$C = 0,3203 \times \Delta DO \quad [\text{g/l of D-Lactic acid in the sample}]$$

$$C = 0,3232 \times \Delta DO \quad [\text{g/l of L-Lactic acid in the sample}]$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 10, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

Assay control

D/L-Lactic acid reagents must be validated with the use of the standard included in the kit. The standard is ready-to-use.

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