R	1 × 60 mL - Phosphate buffer pH 7 Phenol / Amino-4-antipyrine Cholesterol oxidase / Peroxydase
C	1 × 2 mL - Control solution

v028-21031

Sample preparation

The concentration of Cholesterol in the sample used in the assay procedure has to be between 0.05 to 1 g/l.

Precision

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure

Wavelength: 500nm / Optical path: 1cm / Temperature: 20-25 / 37°C
Measurement: against water or air.

	Blank	Sample
Reagent	2,0 ml	2,0 ml
Mix and read DO	DO1	DO1
Water	100 µl	0
Sample	0	100 µl
Mix and read DO after 20 min	DO2	DO2

Calculation

Determine the values for blank and samples:

$$\Delta DO = [DO 2 - DO 1]_{\text{sample}} - [DO 2 - DO 1]_{\text{blank}}$$

Cholesterol concentration is given by:

$$C = \frac{V \times MW}{\varepsilon \times l \times v \times 1000} \times \Delta DO \quad (\text{g/L})$$

In the assay procedure:

$$C = 1,1598 \times \Delta DO \quad [\text{g/l of Cholesterol in the sample}]$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

Assay control

Cholesterol reagents must be validated with the use of the include control.

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