R1 Amidon	1 × 30 mL - Buffer pH 4,5
R2 Amidon	1 × 0,6 mL - Amyloglucosidase
C Amidon	1 × 1 gramme - amidon / starch
Kit #010	1 × 1 kit - D-Glucose

### Sample preparation:

The concentration of Starch in the sample used in the assay procedure has to be between 0.05 to 0.8 g/l after the solubilization and hydrolysis steps.

Starch or starch containing products have to be pretreated before the assay to convert the starch into a soluble form. Because of the simple handling, dimethylsulfoxide DMSO\* is recommended.

### Solubilization of starch with DMSO\* and HCl\*

Weigh 100 mg of sample. Add 20 ml DMSO\* and 5 ml HCl\* 36%. Close flask, and incubate for 30-60 minutes at 60°C with constant stirring. Cool quickly at 20-25°C and add 50 ml water. Adjust to pH 4-5 with sodium hydroxide\* 5M. Transfert in a 100 ml volumetric flask and fill up the mark with water.

### 2 step Assay procedure:

#### 1. Hydrolysis with Starch reagents

Mix together 10 ml R1 and 50 µl R2 (or 20ml/100µl or 30ml/150µl). This preparation, named "Hydrolysis reagent", is stable 4 week if stored at 2-8 °C.

	Reagent blank	Sample blank	Sample
Hydrolysis reagent	1.0 ml	-	1.0 ml
Water	0.5 ml	1.0 ml	-
Sample	-	0.5 ml	0.5 ml

Mix, incubate for 15 min at 60°C

#### 2. Glucose measurement with Glucose kit

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°  
 Measurement: against water or air.

	Reagent blank	Sample blank	Sample
R 1	1.0 ml	1.0 ml	1.0 ml
Water	1.0 ml	1.0 ml	1.0 ml
Hydrolysed sample	0.3 ml	0.3 ml	0.3 ml
Mix and read the DO	DO1	DO1	DO 1
R 2	0.02 ml	0.02 ml	0.02 ml
Mix and read DO after 20 min	DO 2	DO 2	DO 2

**Calculation:** Determine the values

$$\Delta DO = [ DO 2 - DO 1 ]$$

And

$$\Delta \text{Starch} = \Delta DO_{\text{sample}} - \Delta DO_{\text{sample blank}} - \Delta DO_{\text{reagent blank}}$$

$$\text{Starch concentration is given by: } C = \frac{V \times MW}{\epsilon \times l \times v \times 1000} \times \Delta \text{Starch} \quad (\text{g/L})$$

In the assay procedure:

$$C = 0,60462 \times \Delta DO \quad [\text{g/l of starch in the sample sample}]$$

Dilution factor of the sample has to be considered in the calculation.

$$\text{With solid sample: } \text{Starch Content} = \frac{C (\text{g/L})}{\text{weight}} \times 100 \quad (\text{g}/100\text{g})$$

### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

### Assay control

Starch reagents must be validated with the use of the standard included in the kit. The standard is to prepare (solubilisation with DMSO\* and HCl\*).

### Precision:

Under the described conditions, measurement accuracy is 5% on a control solution.

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