R1	1 × 30 mL - Buffer pH 8
R2	30 tablets - NADH tablets
R3	1 × 0,4 mL - NADH Peroxydase
R4	1 × 1,6 mL - Sulfite Oxydase
C	2g - Sulfite Control

v040E-21121

### Reagents preparation:

Dissolve 1 tablet of NADH (R2) with 1 ml of solution R1.

### Sample preparation:

The amount of sulphite present in the reaction cuvette should range between 0.6 and 30 µg of sulphite.

Thus, if a sample volume of 0.10 mL is used the sample solution must be diluted to yield a sulphite concentration between 6 and 300 mg/L.

However, the sample volume can vary from 0.10 to 2.00 mL, by replacing the water. So the analyte range can vary from 0.30 to 300 mg/L depending on the chosen sample volume.

L-Ascorbic acid inhibits sulfite oxidase.

Higher concentrations of L-ascorbic acid in the assay reacts with hydrogen peroxide and thus produce results that are too low : it is necessary to remove L-ascorbic acid with the use of ascorbate oxidase before the sulphite assay.

### Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

### Assay procedure:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°

Measurement: against water or air.

	Blank	Sample
R 1 (+R2)	1.0 ml	1.0 ml
R 3	0.01 ml	0.01 ml
Water	2.0 ml	1.9 ml
Sample	0	0.1 ml
Mix and read the DO	<b>DO1</b>	<b>DO 1</b>
R 4	0.05 ml	0.05 ml
Mix and read DO after 30 min	<b>DO 2</b>	<b>DO 2</b>

### Calculation

Determine the values for blank and samples:

$$\Delta DO = [ DO 1 - DO 2 ]_{\text{sample}} - [ DO 1 - DO 2 ]_{\text{blank}}$$

Sulphite concentration is given by:

$$C = \frac{V_{\text{tot}} \times MW}{\varepsilon \times l \times v_{\text{éch}} \times 1000} \times \Delta DO$$

Depending of the chosen sample volume:

Sample volume	Measurement range	Calcul
0,1 ml (+1,9 ml H <sub>2</sub> O)	6 – 300 mg/L SO <sub>2</sub>	0,31115 x ΔDO
1 ml (+1 ml H <sub>2</sub> O)	0,6 – 30 mg/L SO <sub>2</sub>	0,031115 x ΔDO
2 ml (+0 ml H <sub>2</sub> O)	0,3 – 15 mg/L SO <sub>2</sub>	0,015557 x ΔDO

Dilution factor of the sample has to be considered in the calculation.

### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.

Take necessary precautions for the use of laboratory reagents.

### Assay control

Sulphite kit reagents must be validated with the use of the standard included in the kit. The standard must be prepared :

Dissolve 60 mg (~ 30 mg sulfite) of included powder control in 100 ml of a 1 g/L citric acid solution. Use this 300 mg SO<sub>2</sub> / L solution on the same day of preparation.

Sulphite powder absorb moisture and are oxidized easily and aqueous solutions are instable. Therefore, results of < 100% have to be expected in analysis.

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