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### Preparation of the standard solution:

It is necessary to prepare a standard solution to dilute before use

Use one of the following products and to weigh in a precise way in a graduated flask of 100 ml the quantity indicated

Sodium sulphite, anhydrous, AR grade (Na <sub>2</sub> SO <sub>3</sub> , M=126.04 g/mol, approx 50.8% SO <sub>2</sub> )	60 mg
Sodium bisulphite, AR grade (Na <sub>2</sub> S <sub>2</sub> O <sub>5</sub> , M=190.10 g/mol, approx 67.4% SO <sub>2</sub> )	45 mg
Potassium bisulphite, AR grade (K <sub>2</sub> S <sub>2</sub> O <sub>5</sub> , M=222.33 g/mol, approx 57.6% SO <sub>2</sub> )	50 mg

That corresponds to 300 mg/L of free SO<sub>2</sub>.

Dilute this solution to the 1/2: 0.5 ml + 0.5 ml of demineralised water. The free SO<sub>2</sub> concentration is then **150 mg/L**.

### Preparation of the sample:

The free SO<sub>2</sub> concentration in the sample used for the test has to be between 5 and 250 mg/L.

### Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

### Assay procedure:

Wavelength: 580nm / Optical path: 1cm / Temperature: 20-37°

Measurement: against water or air.

	Blank	Standard	Sample
R 1	0,1 ml	0,1 ml	0,1 ml
Water	1,0 ml	0,985 ml	0,985 ml
Standard	0	0,015 ml	0
Sample	0	0	0,015 ml
Mix and read the DO	<b>DO 1</b>	<b>DO 1</b>	<b>DO 1</b>
R 2	0,1 ml	0,1 ml	0,1 ml
Mix and read the DO after 20 mn	<b>DO 2</b>	<b>DO 2</b>	<b>DO 2</b>

### Calculation

Determine the following values for each cuvette  
**ΔDO = DO 2 – DO 1**

Determine the absorbance differences then

For standard : **ΔA<sub>standard</sub> = ΔDO<sub>standard</sub> – ΔDO<sub>blank</sub>**

For sample : **ΔA<sub>sample</sub> = ΔDO<sub>sample</sub> – ΔDO<sub>blank</sub>**

Determine the free SO<sub>2</sub> concentration in the samples as described below

$$C(\text{mg/L of free SO}_2) = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}}$$

Dilution factor of the sample has to be considered in the calculation.

### Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry;  
Storage conditions : 18 – 25 °C, darkness.

### Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.  
Take necessary precautions for the use of laboratory reagents.

### General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

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