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| | |
|----|-----------|
| RA | 1 × 75 mL |
| RB | 1 × 75 mL |

Preparation of the standard solution:

It is necessary to prepare a standard solution to dilute before use

Use one of the following products and to weigh in a precise way in a graduated flask of 100 ml the quantity indicated

| | |
|--|-------|
| Sodium sulphite, anhydrous, AR grade (Na ₂ SO ₃ , M=126.04 g/mol, approx 50.8% SO ₂) | 60 mg |
| Sodium bisulphite, AR grade (Na ₂ S ₂ O ₅ , M=190.10 g/mol, approx 67.4% SO ₂) | 45 mg |
| Potassium bisulphite, AR grade (K ₂ S ₂ O ₅ , M=222.33 g/mol, approx 57.6% SO ₂) | 50 mg |

That corresponds to 300 mg/L of total SO₂.

Dilute this solution to the 1/6: 1.5 ml + 7.5 ml of demineralised water. The total SO₂ concentration is then 50 mg/L.

Preparation of the sample:

The total SO₂ concentration in the sample used for the test has to be between 5 and 300 mg/L.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 405nm / Optical path: 1cm / Temperature: 20-37°

Measurement: against water or air.

| | Blank | Standard | Sample |
|---------------------------------|----------|----------|----------|
| R 1 | 0,125 ml | 0,125 ml | 0,125 ml |
| Water | 1,250 ml | 1,235 ml | 1,235 ml |
| Standard | 0 | 0,015 ml | 0 |
| Sample | 0 | 0 | 0,015 ml |
| Mix and read the DO | DO 1 | DO 1 | DO 1 |
| R 2 | 0,125 ml | 0,125 ml | 0,125 ml |
| Mix and read the DO after 20 mn | DO 2 | DO 2 | DO 2 |

Calculation

Determine the following values for each cuvette

$$\Delta DO = DO 2 - DO 1$$

Determine the absorbance differences then

$$\text{For standard : } \Delta A_{\text{standard}} = \Delta DO_{\text{standard}} - \Delta DO_{\text{blank}}$$

$$\text{For sample : } \Delta A_{\text{sample}} = \Delta DO_{\text{sample}} - \Delta DO_{\text{blank}}$$

Determine the total SO₂ concentration in the samples as described below

$$C(\text{mg/L of total SO}_2) = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{s tan dard}}} \times C_{\text{s tan dard}}$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry.
Storage conditions: 18 – 25 °C, darkness.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.
Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

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