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
V43-19051

By mixing reduced sample and chromogen in an acidic buffered solution, the reaction produces a measurable blue color within few minutes.

Reagents don't contain acetic acid, so there are no interference for this measurement. In addition, performance has been improved for the analysis of red wines.

Reagents preparation:

Reagent N°1: Reagent R1 + Reagent R3

Mix gently 20 vol. of R1 reagent and 1 vol. of R3 reagent
Ex: 20ml R1 + 1 ml R3

Stability of reagents n°1 : 10 days at +4°C.

Reagent N°2: Reagent R2

The reagent is ready to use.

Stability of reagents: see the expire date of the kit.

Automates use:

Range of measurement: The test was developed to determine the concentration in Iron ranging between 0.1 and 15 mg/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

R1 = 200 / Sample = 40

Incubation 2 min – Lecture 1

R2 = 10

Incubation 3 min – Lecture 2

Multipoint calibration – Linear

End-point

Wavelength: I = 590 nm / II = -

Depending on the used matrix, the calibration could have to be adjusted.

Manual use:

Contact us to get assay protocol.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Assay control

Iron reagents must be validated with the use of a control, include in the kit.

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