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By mixing reduced sample and chromogen in an acidic buffered solution, the reaction produces a measurable blue color within few minutes. Reagents don't contain acetic acid, so there are no interference for this measurement. In addition, performance has been improved for the analysis of red wines.

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Be careful !! New preparation of the working solutions:

Mix 1 mL of R3 reagent in 20 ml of R1 reagent. This solution is stable 10 days at 4°C.

Preparation of the standard solution:

It is necessary to prepare a standard solution, included in the kit

Preparation of the sample:

The iron concentration in the sample used for the test has to be between 0.1 and 15 mg/L.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 590nm / Optical path: 1cm / Temperature: 20-37°
Measurement: against water or air.

	Blank	Standard	Sample
R 1 (with reconstituted R3)	1 ml	1 ml	1 ml
Water	0,2 ml	0	0
Standard	0	0,2 ml	0
Sample	0	0	0,2 ml
Mix and read the DO	DO 1	DO 1	DO 1
R 2	0,05 ml	0,05 ml	0,05 ml
Mix and read the DO after 20 mn	DO 2	DO 2	DO 2

Calculation:

Determine the following values for each cuvette

$$\Delta DO = DO 2 - DO 1$$

Then, determine the absorbance differences

For standard : $\Delta A_{\text{standard}} = \Delta DO_{\text{standard}} - \Delta DO_{\text{blank}}$

For sample : $\Delta A_{\text{sample}} = \Delta DO_{\text{sample}} - \Delta DO_{\text{blank}}$

Finally, determine the Iron concentration in the samples as described below

$$C_{\text{iron}} (\text{mg} / \text{L}) = \left[\frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}} \right] \times F$$

Where F = dilution factor of the sample (if diluted)

Storage instructions and reagent stability:

The reagents are stable up to the indicated date of expiry. Storage conditions: 2 – 8°C, darkness.

Warnings and precautions:

Do not swallow. Avoid contact with the skin and mucous membranes. Take necessary precautions for the use of laboratory reagents.

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