R1	1 x 150 mL - Buffer pH 4.8
R2	10 doses - Reductant
R3	1 x 15 mL - Chromogen
S	1 x 5 mL - Standard solution / Solution standard

V44-19071

The sample is mixed in a pH 4,8 buffered solution containing the chromogen. The reaction develops in a few minutes a colouring measurable red.

The reagents don't contain acetic acid, so there are no interference for this measurement. In addition, no crystals will appear during refrigerated storage or on-board storage.

**Reagents preparation:**

**Reagent N°1: Reagent R1 + Reagent R2**

Mix gently 15 vol. of R1 reagent and 1 dose of R2 reagent  
Ex: 15ml R1 + 1 dose R2

Stability of reagents n°1 : 10 days at +4°C.

**Reagent N°2: Reagent R3**

The reagent is ready to use.

Stability of reagents: see the expire date of the kit.

**Automates use:**

**Range of measurement:** The test was developed to determine the concentration in Copper ranging between 0.25 and 1.5 mg/L in the sample with the use of the automate parameters setting supplied by Biosentec.

**Procedure of test:** For a use on automat, applications are available on request.

Follow the analytical procedure:

- R1 = 150 / Sample = 35
- Incubation 2 min – Lecture 1
- R2 = 15
- Incubation 5 min – Lecture 2
- Multipoint calibration – Linear
- End-point
- Wavelength: I = 580 nm / II = -

Depending on the used matrix, the calibration could have to be adjusted.

**Manual use:**

Contact us to get assay protocol.

**Storage instructions and reagent stability**

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

**Warnings and precautions**

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

**Assay control**

Copper reagents must be validated with the use of a control, include in the kit.

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