R1	1 × 150 mL - Buffer pH 4.8
R2	10 doses - Reductant
R3	1 × 15 mL - Chromogen
S	1 × 5 mL - Standard solution / Solution standard

The sample is mixed in a pH 4,8 buffered solution containing the chromogen. The reaction develops in a few minutes a colouring measurable red.

The reagents don't contain acetic acid, so there are no interference for this measurement. In addition, no crystals will appear during refrigerated storage or on-board storage.

Preparation of the standard solution:

It is necessary to prepare a standard solution, include in the kit

Preparation of the sample:

The copper concentration in the sample used for the test has to be between 0.25 and 1.5 mg/L.

Precision:

Under the described conditions, measurement accuracy is 5% on a control solution

Assay procedure:

Wavelength: 580nm / Optical path: 1cm / Temperature: 20-37°

Measurement: against water or air.

	Blank	Standard	Sample
R 1	1 ml	1 ml	1 ml
Water	0,2 ml	0	0
Standard	0	0,2 ml	0
Sample	0	0	0,2 ml
Mix and read the DO	DO 1	DO 1	DO 1
R 2	0,1 ml	0,1 ml	0,1 ml
Mix and read the DO after 20 mn	DO 2	DO 2	DO 2

Calculation

Determine the following values for each cuvette

$$\Delta DO = DO 2 - DO 1$$

Determine the absorbance differences then

For standard : $\Delta A_{\text{standard}} = \Delta DO_{\text{standard}} - \Delta DO_{\text{blank}}$

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determine the Copper concentration in the samples as described below

$$C(\text{mg/L of Copper}) = \frac{\Delta A_{\text{sample}}}{\Delta A_{\text{standard}}} \times C_{\text{standard}}$$

Dilution factor of the sample has to be considered in the calculation.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry.

Storage conditions: 18 – 25 °C, darkness.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.

Take necessary precautions for the use of laboratory reagents.

General information on sample preparation

- Use liquid samples, clear, colorless and with pH practically 8, directly, or after dilution.
- Crush or homogenize solid or semi-solid samples
- Deproteinize samples containing proteins with perchloric acid; alternatively clarify with Carrez reagents.
- Extract samples containing fat with hot water

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