Amino Nitrogen + NAC + OPA \longrightarrow Isoindole derivative

RA	1 \times 100 mL – Buffer
RB	1 \times 100 mL – NAC / OPA
C	1 \times 1 gr – Iso-Leucine

v48-19191

In basic medium, primary amino groups are derivatized by OPA and NAC to form a derivative, with a maximal absorbance at 340 nm. The optical density measured at 340nm is directly proportional to the Primary amino nitrogen in the sample.

Reagents preparation:

Reagent N°1: Reagent A

Reagent N°2: Reagent B

The reagents are ready to use.

Stability of reagents: see the expire date of the kit.

Automates use:

Range of measurement: The test was developed to determine the concentration in Primary amino nitrogen ranging between 20 and 300 mg/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

R1 = 135 / Sample = 4

Incubation 2 min – Lecture 1

R2 = 135

Incubation 5 min – Lecture 2

Multipoint calibration – Linear

End-point

Wavelength: I = 340 nm / II = 700-800 nm

Depending on the used matrix, the calibration could have to be adjusted.

Manual use:

Contact us to get assay protocol.

Storage instructions and reagent stability

Reagents are light sensible.

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Standard curve and Control preparation :

Prepare one control at 100 mg/L of primary amino nitrogen.

Weight 93.7 mg of isoleucine, and dissolve it in 100 ml water.

Prepare the standard curve with a solution at 300 mg/L of primary amino nitrogen.

Weight 281 mg of isoleucine, and dissolve it in 100 ml water.

Solutions are stables 15 days, and must to be store at 2-8°C.

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