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In basic medium, primary amino groups are derivatized by OPA and NAC to form a derivative, with a maximal absorbance at 340 nm. The optical density measured at 340nm is directly proportional to the Primary amino nitrogen in the sample.

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Reagents preparation:

Reagents are ready to use.

Reagents must be store between 2 and 8°C.
Reagents are light sensible.

Standard curve and Control preparation :

Prepare one control at 100 mg/L of primary amino nitrogen. Weight 93.7 mg of isoleucine, and dissolve it in 100 ml water.

Prepare the standard curve with a solution at 300 mg/L of primary amino nitrogen. Weight 281 mg of isoleucine, and dissolve it in 100 ml water.

Solutions are stables 15 days, and must to be store at 2-8°C.

Automates use:

Range of measurement: The test was developed to determine the primary amino nitrogen ranging between 20 and 300 mg/L of.

Precision: Under the described conditions, measurement accuracy is 5% on a control solution

Procedure of test: For a use on automat, programmings are available on request.

Assay procedure for a manual use:

Wavelength: 340nm / Optical path: 1cm / Temperature: 20-37°
Measurement: against water or air.

	Blank	Standard	Sample
R 1	0,75 ml	0,75 ml	0,75 ml
R 2	0,75 ml	0,75 ml	0,75 ml
Mix and read DO		DO 1	DO 1
Water	0,02 ml	–	–
Standard	–	0,02 ml	–
Sample	–	–	0,02 ml
Mix and read DO at 2 min	DO 2	DO 2	DO 2

Calculation

Determine the values for samples:

$$\Delta DO = [DO 2 - DO 1]_{\text{sample}} - [DO 2 - DO 1]_{\text{blank}}$$

Trace the standard curve and report your sample.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored between 2 and 8°C.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.

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