Ca⁺⁺ + Arsenazo III → Blue complex

R1	1 × 100 mL	Buffer pH 7 / 8-hydroxyquinoline-5-sulfonic acid / Arsenazo III
S	1 × 5 mL	Standard solution / Solution standard

V50-19071

By using 8-hydroxyquinoline-5-sulfonic acid to eliminate the interference of magnesium, calcium ions combine with Arsenazo III to produce a blue colored complex at a neutral solution. The absorbency increase is directly proportional to the concentration of calcium.

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Reagents preparation:

The reagent is ready to use.
Stability of reagents: see the expire date of the kit.

Automates use:

Range of measurement: The test was developed to determine the concentration in Calcium ranging between 20 and 200 mg/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

R1 = 200
Incubation 2 min – Lecture 1
Sample = 2
Incubation 3 min – Lecture 2
Multipoint calibration – Linear
End-point

Wavelength: I = 660 nm / II = 700 nm

Depending on the used matrix, the calibration could have to be adjusted.

Manual use:

Contact us to get assay protocol.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Assay control

Calcium reagents must be validated with the use of a control.

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