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The reagent is ready to use.
Stability of reagents: see the expire date of the kit.

Automates use:

Range of measurement: The test was developed to determine the concentration in Glucose / Fructose ranging between 0.1 and 8 g/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

- R1 = 250 / Sample = 3
- Incubation 2 min – Lecture 1
- R2 = 50
- Incubation 5 min – Lecture 2
- Multipoint calibration – Linear
- End-point
- Wavelength: I = 340 nm / II = 700 nm

Depending on the used matrix, the calibration could have to be adjusted.

Manual use:

Contact us to get assay protocol.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.
Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Standard and assay control

Glucose / Fructose reagents can be calibrated and validated with the use of a synthetic solution, available on catalogue :

210 : Acid Sugars multi-control (GF 8.0 g/L) – 25ml

Reference :

Resolution OIV-OENO 600-2018

For a better stability, TEA buffer was replaced by TRIS buffer.

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