RA	1 × 100 mL – Buffer pH8.4 / Tampon pH8,4
RB	4 x ATP/CoA/NAD – To reconstitute 4 x ATP/CoA/NAD – A reconstituer
RC	1 × 1 mL – Enzymes
RD	1 × 20 ml – Dilution solution / Solution de dilution

Reagents preparation:

Reagent N°1: 25 ml Reagent A + 1 tube of Reagent B

Take 25 ml of Reagent A.

Dissolve the reagent B with 3 ml of this reagent A. Mix.

Entirely pour the reagent thus reconstituted in the remainder of reagent A taken beforehand.
Mix.

Stability of reagent N°1: 10 days.

Reagent N°2: 0,25 ml Reagent C + 5 ml Reagent D

Swirl bottles to mix contents before use.

Take 0,25 ml Reagent C and 5 ml Reagent D.

Gently mix by inversion (10 times) these 2 reagents.
Let rest the solution 5 minutes before use.

Stability of reagent N°2: 10 days

Automates use:

Range of measurement: The test was developed to determine the concentration in acetic acid ranging between 0.1 and 1 g/L in the sample with automate parameters setting supplied by Biosentec.

Precision: Under the described conditions, measurement accuracy is 5% on a control solution

Procedure of test: For a use on automat, programmings are available on request.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.

Assay control

Acetic acid reagents must be validated with the use of the control, available on catalogue.

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