RA	1 × 100 mL – TRIS pH7.6 / PEP / MgCl ₂ / KCl / PVP
RB	1 × 5 mL – ATP / NADH / Stabilizers
RC	1 × 1 mL – L-Lactate dehydrogenase / Stabilizers
RD	1 × 25 mL – Acetate Kinase / Pyruvate Kinase / Stabilizers

v77-19511

Reagents preparation:

Reagent N°1: Reagent A + Reagent B + Reagent C

Mix gently 20 vol. of RA reagent + 1 vol. of RB reagent + 0.2 vol. of RC reagent

Wait for 15 minutes before use.

Stability of reagents n°1 : 10 days.

Reagent N°2: Reagent D

Reagent D is ready to use.

Stability of reagents n°2 : see expiration date of reagent RD.

Automates use:

Range of measurement: The test was developed to determine the concentration in Acetic acid ranging between 0.05 and 1 g/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

R1 = 250 / Sample = 3

Incubation 2 min – Lecture 1

R2 = 50

Incubation 5 min – Lecture 2

Multipoint calibration – Linear

End-point

Wavelength: I = 340 nm / II = 700 nm

Depending on the used matrix, the calibration could have to be adjusted.

Manual use:

Contact us to get assay protocol.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Assay control

Acetic acid reagents must be validated with the use of a control, available on catalogue.

Reference :

Resolution OIV-OENO 621-2019

For a better stability, MOPS buffer was replaced by TRIS buffer.

Reagents formulation eliminates the influence of pyruvate in acetate results.

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