Liquid reagent / Réactif liquide	
RA	1 × 100 mL – Buffer pH10 / Glutamate
RB1	1 × 30 mL – NAD / Stabilizer
RB2	1 × 1,5 mL – GPT / D-LdH

v079-23031

Reagents preparation:

Reagent N°1: Reagent 1
Reagent is ready to use

Stability of reagents n°1 : see the expire date.

Reagent N°2: Reagent B1 + Reagent B2
Mix gently 2 vol. of RB1 and 0.1 vol. of RB2
Wait for 15 minutes before use.

Stability of reagents n°2 : 30 days.

Automates use:

Range of measurement: The test was developed to determine the concentration in D-lactic acid ranging between 0.020 and 1.000 g/L in the sample with the use of the automate parameters setting supplied by Biosentec.

Procedure of test: For a use on automat, applications are available on request.

Follow the analytical procedure:

- R1 = 200 / Sample = 2
- Incubation 2 min – Lecture 1
- R2 = 40
- Incubation 5 min – Lecture 2
- Multipoint calibration – Linear
- End-point
- Wavelength: I = 340 nm / II = 700 nm

Depending on the used matrix, the calibration could have to be adjusted.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Onboard stability: 5 days

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes.

See the material security data sheet for more information.

Standard and Assay control

D-lactic acid reagents can be calibrated and validated with the use of a synthetic solution, available on catalogue :

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