R1A	1 x 24 mL Buffer pH6.5 / HexoKinase / PhosphoGlucosomérase / PVP / stabilisants
R1B	1 x 0.3 mL Maltase / Stabilisants
R2A	1 x 30 mL Buffer pH7.5 / MgSO ₄ / ATP / NADP
R2B	1 x 6 mL Hexokinase / Phosphoglucoisomérase / Glucose-6-phosphate dehydrogenase / Stabilisants

v70-23031

Reagent's preparation:

Reagent N°1: Reagent R1A + Reagent R1B

Mix gently 8 ml of R1A reagent and 0.1 ml reagent R1B
 Stability of reagents: 30 days.

Reagent N°2: Reagent R2A + Reagent R2B

Mix gently 10 ml of R2A reagent and 2 ml reagent R2B
 Stability of reagents: 30 days

Measuring range:

The test was developed to determine the concentration of fermentable sugars: Maltose / Maltotriose / Sucrose / Glucose / Fructose ranging between 0.1 and 10 g/L in the sample.

Automates use:

Applications are available on request.

Follow the analytical procedure:

R1 = 120 / Sample = 2

Incubation 5 min, 37°C – Lecture 1

R2 = 180

Incubation 5 min ; 37 °C – Lecture 2

Multipoint calibration – Linear

Endpoint

Wavelength: I = 340 nm / II = 700 nm

Depending on the used matrix, the calibration could have to be adjusted. Contact us to get assay protocol for manual use.

Specificity:

α -glucosidase only hydrolyzes α -glucosidic bonds. Lactose is therefore not hydrolyzed.

Maltose, Sucrose and Maltotriose are hydrolysed. Under the test conditions, the maltotetraose hydrolysis rate is less than 5%; Maltopentaose, larger dextrans and starch are not hydrolyzed. The determination of D-Glucose is specific.

Storage instructions and reagent stability

The reagents are stable up to the indicated month of expiry, if stored at 2-8 °C, contamination is avoided.

Warnings and precautions

Do not swallow. Avoid contact with the skin and mucous membranes. See the material security data sheet for more information.

Standard and assay control

Fermentable sugars reagents can be calibrated and validated with the use of a synthetic solution, available on catalogue :

226 : Fermentable sugars (maltose+sucrose = 10.0 g/L equivalent glucose) – 25ml

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